


## VieLight Neuro

This device stimulates your brain safely and non-invasively with pulsed near-infrared light.

### How does VieLight Neuro work?

Powerful, yet simple, the VieLight is a photobiomodulation (aka cold laser therapy) head cap, utilizing multiple frequencies of red light to help enhance a myriad cellular functions. The transcranial LED diodes are designed to target the default mode network directly, and the intranasal laser stimulates your brain's prefrontal cortex (decision-making) and olfactory bulb (memory) for optimal effectiveness.

The pulse rate of this device can be set to Alpha or Gamma format. Alpha induces a calming effect with improved recovery of neuroinflammation, while Gamma enhances memory, cognition, awareness, and mental focus.



The nasal passage is a direct gateway for light energy to reach your brain's prefrontal cortex and olfactory bulb. This unique approach sets the VieLight apart from other devices in its class.



### Clinically Proven Benefits

- Increases synaptogenesis
- Increases angiogenesis
- Increases blood flow
- Anti-inflammatory
- Increases neuron progenitor cells
- Increases antioxidants
- Increases DNA protection factors
- Increases neurotrophins
- Reduces neuronal excitotoxicity



### Symptoms Addressed

- Concussion/PCS
- Fatigue
- Headaches
- Depression
- Anxiety
- Memory challenges
- Dyssynergia
- Brain fog